

Dr. Q 分子生物試劑

Dr.Q Taq DNA Polymerase

Product name: Dr.Q Taq DNA

Polymerase,500u

Catalogue no.:Dr.QPCR004

Cat No	Package	conc
Dr.QPCR004	500u	5U/ul

Description:

DR.Q Taq DNA polymerase is a thermostable enzyme of isolated from *Thermus aquaticus*. This enzyme contain 5'-3' polymerase and 5'-3' exonuclease activity.

Storage buffer:

50mM Tris-HCl pH7.9, 50mM KCl, 0.1mM EDTA, 1mM DTT, 0.5mM PMSF , 50% glycerol.

10X reaction buffer: buffer A containing 15mMMgCl₂.

buffer B without MgCl₂

Unit description:

one unit is defined as the amount of enzyme that will incorporate 10n mole of dNTP into acid-insoluble material in 30 minutes at 74oC. The reaction conditions are: 50mM Tris-HCl pH8.8, 50mM NaCl, 5mM MgCl₂, 200uM each of dATP, dCTP, dGTP, dTTP, 10 ug activated calf thymus DNA and 0.1mg/ml BSA in a final volume of 50 ul.

Storage:

50% glycerol (v/v), 20 mM Tris-HCl pH 8.7 at -20°C, 100 mM KCl, 0.1 mM EDTA.

Source: E coli clone

Quality control:

The enzyme is free of nicking and priming activities, exonucleases and non-specific endonucleases. SDS/PAGE - 95 kD band. Activity and stability tested via thermo-cycling. The error rate per nucleotide per cycle is $\sim 2.5 \times 10^{-5}$; the accuracy is $\sim 4 \times 10^4$. Estimated half life at 95°C is 1.5 hours.

Shipping and Storage conditions:

Shipping and temporary storage at -20 for up to 1 month at room temperature has no detrimental effects on the quality of DR.Q Taq DNA polymerase.

尉徠生物科技有限公司

FUTURE SCIENTIFIC CO.,LTD

地址：桃園縣龜山鄉文化一路259號工學大樓1F創新育成中心5884室

TEL：(03)211-8800#5884；(03)211-8299 MOBILE：0955-103293

Dr. Q 分子生物試劑

PCR reaction mix:

Component	Volume
DR.Q Taq	0.5-1ul
10X buffer	10 ul
10mM dNTP	2 ul
Primer1 (20 pmol)	2-4 ul
Primer2 (20 pmol)	2-4 ul
template	1-10 ul
ddH ₂ O	Up to 100 ul
Total	100 ul

PCR cycles

Step	Temperature	Time	Cycle
Initial denaturation	94-95°C	1-3 mins	1
Denaturation	94-95°C	10-60sec	25-35
Annealing	50-68°C	10-30sec	
Extension	72°C	1min/1kb	
Final extension	72°C	1-10 mins	1

IMPORTANT:

Annealing temperature should be 2-6°C lower than the primer melting temperature.

尉徠生物科技有限公司

FUTURE SCIENTIFIC CO.,LTD

地址：桃園縣龜山鄉文化一路259號工學大樓1F創新育成中心5884室

TEL：(03)211-8800#5884；(03)211-8299 MOBILE：0955-103293