

## 100 bp DNA Ladder

Cat No. FM100B  
Pack size: 500  $\mu$ l  
Conc: 500  $\mu$ g

### Description

100 bp DNA Ladder is suitable for size of linear double stranded DNA fragments from 100 bp to 3 kb. The 13 bands of the ladder contain fragments following sizes: 100, 200, 300, 400, 500, 600, 700, 800, 900, 1000, 1500, 2000 and 3000 bp. The intensity of 500 bp is about three times than the other bands to serve as reference points.

### Storage buffer

10 mM Tris-HCl (pH8.0), 1 mM EDTA

### Storage and shipping

For long term storage store at  $-20^{\circ}\text{C}$   
DNA ladder is stable for at least 6 months at  $4^{\circ}\text{C}$

## Gel analysis

