



## BioDiamond dNTP mix

Cat. No. DMDD010

DMDD011.

DMDD012

Cat No	Products	Package
DMDD010	2.5 mM dNTP mix	1ml
DMDD011	10 mM dNTP mix	1ml
DMDD012	100 mM dNTP mix	1mlx4

### Description

The deoxynucleotides (PCR grade) are suitable for many applications where high quality reagents are required. Such procedures include RT, PCR, RT-PCR, DNA labeling reaction and sequencing/ cycle Sequencing analysis.

### Purity assay

HPLC analysis >99%

NMR analysis (inorganic phosphates)–passed Exo-endo deoxyribonucleases test–passed UV analysis–passed Spectrophotometry–passed.

### Functional assays

amplification of 8 kb PCR fragment from genomic DNA with Taq DNA polymerase–passed amplification of 0.6 kb PCR fragment from genomic DNA with Pfu DNA Polymerase -passed

### Storage

store at –20°C

蔚徠生物科技有限公司  
FUTURE SCIENTIFIC CO.,LTD  
地址：桃園市龜山區文化一路259號工學大樓1F創新育成中心5884室  
TEL：(03)211-8800#5884；(03)211-8299 MOBILE：0955-103293