



BioDiamond RNA Reagent

Cat. No. DMDRT006

Pack size: 200ml

Description

RNA REAGENT is a ready-to-use reagent for the isolation of total RNA from cells and tissues. Just a single-step is used during cell homogenization or lysate. 1ml of reagent will treat small quantities of tissue (50-100 mg) or cells (5×10^6) well.

Storage condition

2 to 8°C stable for 12 months


Application

RT-PCR

Northern hybridization

RNase protection d. poly-A+ RNA selection Differential display and microarray assay

Quality result



a variety of RNA species has an A260/A280 ratio 1.8 to 2.1 when diluted into TE.

蔚徠生物科技有限公司
FUTURE SCIENTIFIC CO.,LTD
地址：桃園縣龜山鄉文化一路259號工學大樓1F創新育成中心5884室
TEL：(03)211-8800#5884；(03)211-8299 MOBILE：0955-103293